

RECKNEY Instructions for Use

Product: Nasal Strip Refill Tabs

Quantity: 120 pcs

Intended use

Adhesive nasal strip refill tabs intended for external application on the nose to support proper positioning of compatible nasal strips. For external use only.

Directions for use

1. Clean and dry the nose area.
2. Apply one nasal strip refill tab to the skin in the desired position.
3. Ensure it is firmly attached before use.
4. Remove gently after use.

Warnings

Do not use on broken, irritated or infected skin. Discontinue use if irritation or discomfort occurs. Keep out of reach of children. Single-use product. Do not reuse.

Storage

Store in a cool, dry place away from direct sunlight.

Manufacturer

RECKNEY

support@reckney.co.uk